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● References

- [1] Yoshida M, Kadokami T, Momii H, Hayashi A, Urashi T, Narita S, et al. (2012) Enhancement of cardiac performance by bilevel positive airway pressure ventilation in heart failure. *Journal of Cardiac Failure* 18 (12): pp. 912–918.
- [2] 急性心不全治療ガイドライン(2011年改訂版)
- [3] Parati G, Lombardi C, Castagna F, Mattaliano P, Filardi PP, Agostoni P, et al (2016) Heart failure and sleep disorders. *Nature Reviews Cardiology* 13: pp. 389–403.
- [4] Kasai K, Bradley TD (2011) Obstructive Sleep Apnea and Heart Failure. *Journal of the American College of Cardiology* 57 (2): pp. 119–127.
- [5] Tamura A, Kawano Y, Naono S, Kotoku M, Kadota J. (2007) Relationship between beta-blocker treatment and the severity of central sleep apnea in chronic heart failure. *Chest* 131 (1): pp. 130–135.
- [6] Lamba J, Simpson CS, Redfearn DP, Michael KA, Fitzpatrick M, Baranchuk A (2011) Cardiac resynchronization therapy for the treatment of sleep apnoea: a meta-analysis. *Europace* 13 (8): pp. 1174–1179.
- [7] Mortara A, Sleight P, Pinna GD, Maestri R, Capomolla S, Febo O, et al. (1999) Association between hemodynamic impairment and Cheyne-Stokes respiration and periodic breathing in chronic stable congestive heart failure secondary to ischemic or idiopathic dilated cardiomyopathy. *American Journal of Cardiology* 84 (8): 900–904.
- [8] Sasayama S, Izumi T, Matsuzaki M, Matsumori A, Asanoi H, Momomura S, et al. (2009) Improvement of quality of life with nocturnal oxygen therapy in heart failure patients with central sleep apnea. *Circulation Journal* 73 (7): pp. 1255–1262.
- [9] Nakao YM, Ueshima K, Yasuno S, Sasayama S (2016) Effects of nocturnal oxygen therapy in patients with chronic heart failure and central sleep apnea: CHF-HOT study. *Heart and Vessels* 31 (2): pp. 165–172.
- [10] Yoshihisa A, Suzuki S, Miyata M, Yamaki T, Sugimoto K, Kuni H, et al. (2012) 'A single night' beneficial effects of adaptive servo-ventilation on cardiac overload, sympathetic nervous activity, and myocardial damage in patients with chronic heart failure and sleep-disordered breathing. *Circulation Journal* 76 (9): pp. 2153–2158.
- [11] Yoshihisa A, Suzuki S, Yamaki T, et al. Impact of adaptive servo-ventilation on cardiovascular function and prognosis in heart failure patients with preserved left ventricular ejection fraction and sleep-disordered breathing. *Eur J Heart Fail.* 2013;15(5):543-550.
- [12] Yoshihisa A, Suzuki S, Yamauchi H, et al. Beneficial Effects of Positive Airway Pressure Therapy for Sleep-Disordered Breathing in Heart Failure Patients With Preserved Left Ventricular Ejection Fraction. *Clin Cardiol.* 2015;38(7): 413-421.

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● References

- [11] Takama N, Kurabayashi M (2013) Safety and efficacy of adaptive servo-ventilation in patients with severe systolic heart failure. *Journal of Cardiology* 63 (4): pp. 302–307.
- [12] Takama N1, Kurabayashi M (2012) Effect of adaptive servo-ventilation on 1-year prognosis in heart failure patients. *Circulation Journal* 76 (3): pp. 661–667.
- [13] Cowie MR, Woehrle H, Wegscheider K, Angermann C, d'Ortho MP, Erdmann E, et al. (2015) Adaptive Servo-Ventilation for Central Sleep Apnea in Systolic Heart Failure. *The New England Journal of Medicine* 373 (12): pp. 1095–1105.

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